

USER GUIDE

ConnectLine Office Phone Kit



Introduction

The use of office phones with hearing instrument solutions can be difficult, because different brands of office phones use different proprietary/digital encoded signals. Consequently, a traditional analog phone interface - such as the ConnectLine Phone Adapter - will not work like it does with domestic/landline phones.



The only uniform way to interface to these systems across different brands/vendors is through an analog headset interface which either works in parallel with the handset on the phone or, if available, a dedicated headset interface/headset port.

The ConnectLine Office Phone Kit, which includes a Sennheiser Communications UI770 interface box and a special ConnectLine interface cable, enables ConnectLine users to use PABX/digital office phones in two different ways:

1. Wired connection directly to the Streamer.
2. Wireless connection using the ConnectLine Phone Adapter 2.0.



UI770 Interface Box



ConnectLine Interface Cable



Streamer connected directly to the UI770/Office phone

Using the UI box directly connected to the Streamer


The Streamer can be connected directly to the UI770 box, so that ConnectLine users (Streamer Pro 1.0 or later) may benefit from “plug-and-play” headset functionality.

1. Connect the UI 770 box to the office phone as described in the UI770 manual (A).
2. Connect the Streamer directly to the UI770 using the included ConnectLine interface cable (B). The Mini Jack is connected to the Streamer Mini Jack in/out, and the RJ11 connector to the headset socket on the front of the UI770.
3. Connection to the UI770 is established by pressing the AUX selector on the Streamer (C). All handling of the call is done through the UI770 box (or the office phone). Please refer to the UI770 instructions.

Please note!

- ▶ No call alerts will be heard in the Hearing Instruments (since it is a pure analog interface to the office phone)
- ▶ If a cell phone (or ConnectLine Phone Adapter) is paired to the Streamer, incoming calls to either will have priority over the UI770 connection

For more information please visit www.oticonusa.com/connectline



<----->
up to 30 meters
(maximum with free line of sight)

Wireless using Phone Adapter 2.0

Using ConnectLine Phone Adapter 2.0 with the UI770 for a wireless headset connection

Oticon ConnectLine users can benefit from wireless headset functionality by connecting the ConnectLine Phone Adapter 2.0 to the UI770 box.

- 1.** Connect the UI770 box to the office phone as described in the UI770 manual (A).
- 2.** Connect the Phone Adapter 2.0: ConnectLine socket (B) to the UI770 headset port (C) using the included ConnectLine interface cable.
- 3.** Connection to the UI box/office phone is established by pressing the phone button on the Streamer (for pairing details, please refer to the Instructions For Use included with the Phone Adapter). All handling of the call is done through the UI770 box (or the office phone). Please refer to the UI770 instructions.

Please note!

- ▶ No call alerts will be heard in the Hearing Instruments because it is a pure analog interface between the Phone Adapter 2.0 and the office phone

For more information please visit www.oticonusa.com/connectline

People First

People First is our promise
to empower people
to communicate freely,
interact naturally and
participate actively



00001286608000001

oticon
PEOPLE FIRST