



Technical data sheet



	Oticon Opn 1	Oticon Opn 2	Oticon Opn 3	
Speech Understanding	OpenSound Navigator™	Level 1	Level 2	Level 3
	Balancing power effect	100%	50%	50%
	Max. noise removal	9 dB	5 dB	3 dB
	Speech Guard™ LX	Level 1	Level 2	Level 3
	Spatial Sound™ LX	4 estimators	2 estimators	2 estimators
	Soft Speech Booster LX	•	•	•
	Speech Rescue™ LX	•	•	•
Sound Quality	Clear Dynamics	•	•	-
	Spatial Noise Management	•	•	-
	Fitting Bandwidth*	10 KHz	8 KHz	8 KHz
	Processing Channels	64	48	48
	Bass Boost (streaming)	•	•	•
Listening Comfort	Transient Noise Management	4 configurations	On/Off	On/Off
	Feedback shield LX	•	•	•
	Wind Noise Management	•	•	•
	Binaural Coordination***	•	•	•
Personalization & Optimizing Fitting	YouMatic™ LX	3 configurations	2 configurations	1 configuration
	Fitting Bands	16	14	12
	Multiple Directionality Options	•	•	•
	Adaptation Management	•	•	•
	Oticon Firmware Updater	•	•	•
	Fitting Formulas	VAC+, NAL-NL1+2, DSL v5.0	VAC+, NAL-NL1+2, DSL v5.0	VAC+, NAL-NL1+2, DSL v5.0
	Acoustic Notifications	•	•	•
Connecting to the World	Stereo streaming (2.4 GHz)	○	○	○
	Oticon ON App	○	○	○
	ConnectClip	○	○	○
	Remote Control 3.0	○	○	○
	TV Adapter 3.0	○	○	○
	Autophone	○	○	○
	Tinnitus SoundSupport™***	•	•	•
Battery life, hours**	55-60 / 105-115	55-60 / 105-115	55-60 / 105-115	

* Bandwidth accessible for gain adjustments during fitting
 ** Battery size 312 - IEC PR41 / Battery size 13 - IEC PR48.
 Real usage battery life is shown as an estimated interval based on mixed use cases with variable amplification settings and variable input levels, incl. direct stereo streaming from a TV (25% of the time) and streaming from a mobile phone (6% of the time).
 *** If push button is choosen
 • Default ○ Optional - Not included

OTICON | Opn ITC, ITE HS & FS 90



Oticon Opn™ ITC, ITE HS & FS introduce an updated faceplate design.

OpenSound Navigator™ provides better speech understanding by continuously analyzing the environment, balancing all sound sources and attenuating the dominating noise.

TwinLink™ wireless technology combines binaural communication and 2.4 GHz connectivity in stereo directly to external digital devices with very low power consumption. 2.4 GHz is an optional.

Oticon Opn is a Made for iPhone® hearing aid.

Oticon Opn is built on the Velox™ platform, providing frequency resolution in 64 channels (Opn 1).

Fully programmable with updatable firmware, the Velox platform is ready for the future.



For information on compatibility, please visit www.oticon.com




Technical data Measured according to		Ear Simulator IEC 60118-0:1983/AMD1:1994, IEC 60118-0:2015, IEC 60118-1:1995+AMD1:1998 CSV and IEC 60318-4:2010			ZCC Coupler ANSI S3.22-2014, IEC 60118-0:2015 and IEC 60318-5:2006			
Oticon Opn ITC ITE HS & FS 90		Opn 1	Opn 2	Opn 3	Opn 1	Opn 2	Opn 3	
Frequency range Hz		110-9500	110-7500	110-7500	100-7900	100-7500	100-7500	
OSPL90	Peak	131 dB SPL			120 dB SPL			
	1600 Hz	123 dB SPL			115 dB SPL			
	HFA-OSPL90	124 dB SPL			116 dB SPL			
Full-on gain*	Peak	65 dB			55 dB			
	1600 Hz	57 dB			48 dB			
	HFA-FOG	58 dB			50 dB			
Reference test gain		48 dB			39 dB			
Telecoil output (1600 Hz)	1 mA/m field	87 dB SPL			-			
	10 mA/m field	107 dB SPL			-			
	SPLITS L/R	-			96/96 dB SPL			
Total harmonic distortion (Input 70 dB SPL)	500 Hz	2 %			< 2 %			
	800 Hz	2 %			< 2 %			
	1600 Hz	2 %			< 2 %			
Equivalent input noise level	Omni	18 dB SPL			15 dB SPL			
	Dir	28 dB SPL			27 dB SPL			
Battery consumption**	Typical	1.8 mA			1.8 mA			
	Quiescent	1.7 mA			1.7 mA			
Battery life, calculated, hours 312 and 13***		100 / 175			100 / 170			
IRIL (IEC 60118-13:2016)		700/1400/2000 MHz: 20/12/6 dB SPL						

* Measured with the gain control of the hearing aid set to its full-on position minus 20 dB and with an input SPL of 70 dB. This is to obtain a gain response equal to the full-on gain response from e.g. IEC 60118-0+A1:1994 but without influence of feedback.

** Battery current is measured according to IEC 60118-0:1983/AMD1:1994 §7.11, IEC 60118-0:2015 §7.7 and ANSI S3.22:2014 §6.13 after a settling time of minimum 3 minutes.

*** Based on the standardised battery consumption measurement (IEC 60118-0:1983/AMD1:1994). The actual battery life depends on battery quality, use pattern, active feature set, hearing loss and sound environment.



Operating conditions Temperature: +1°C to +40°C Relative humidity: 5% to 93%, non-condensing	Storage and transportation conditions Temperature and humidity should not exceed the following limits for extended periods during transportation and storage. Temperature: -25°C to +60°C Relative humidity: 5% to 93%, non-condensing
--	--