

New Features and Functionality

Opn 1	Premium
Velox™ Platform	•
OpenSound Navigator™	•
TwinLink™ with NFMI and 2.4 GHz	•
Spatial Sound™ LX	•
YouMatic™ LX	•
Speech Guard™ LX	•
Clear Dynamics	•
Wind Noise Management	•
Software	Genie 2
Wireless programming	Fitting Link 3
Cable Fitting	Flex Connect and Cable #3
Made for iPhone®	•
TV Adapter 3.0	•
Remote Control 3.0	•
Battery Size	312
Battery life**	60-65 /55-65

** Real usage battery life is shown as an estimated interval based on mixed use cases with variable amplification settings and variable input levels, incl. direct stereo streaming from a TV (25% of the time) and streaming from a mobile phone (6% of the time).

A new, more discreet miniRITE


Oticon Opn™ miniRITE has a new, discreet design with a smart single push button for easy operation of volume and programs.

Oticon Opn miniRITE offers clients a discreet hearing aid with new features and functionalities including 2.4 GHz wireless technology and Made for iPhone® capabilities.*



miniFit receivers

Select between two receiver types with different output performance, in lengths 0-5



Earpieces

miniFit earpieces include standard domes, MicroMold, and LiteTip. Custom earpieces are available in acrylic, soft silicone or VarioTherm material.



*Oticon Opn is compatible with iPhone SE, iPhone 6s Plus, iPhone 6s, iPhone 6 Plus, iPhone 6, iPhone 5s, iPhone 5c, iPhone 5, 9.7-inch iPad Pro, 12.9-inch iPad Pro, iPad Air 2, iPad Air, iPad (4th generation), iPad mini 4, iPad mini 3, iPad mini 2, iPad mini, and iPod touch (5th and 6th generation). Devices must be running iOS 9.3 or later.

IP68

Proven IP68 water and dust resistant hearing aids



oticon
PEOPLE FIRST