

**Power
your brain**



**New Power
Solutions**



**Alta2 | Nera2 | Ria2
Oticon | Dynamo**

oticon
PEOPLE FIRST

Inium Sense Performance Line Matrix



Features	Alta2 Pro	Alta2	Nera2 Pro	Nera2	Ria2 Pro	Ria2
Custom Styles: IIC through Full Shell*	Yes (except non-wireless IIC & CIC)	Yes (except IIC)	Yes	Yes (except IIC)	Yes	Yes (except IIC)
New BTE Plus Power Design BTE Styles: designRITE, RITE, miniRITE, BTE 13, BTE 13P, BTE 13 Plus Power	Yes	Yes (except designRITE)	Yes	Yes (except designRITE)	Yes	Yes (except designRITE)
Binaural Processing (compression)	Yes (Premium)	No	Yes (Advanced)	No	No	No
Binaural Synchronization (automatics)	Yes	Yes	Yes	Yes	Yes	No
Binaural Coordination (VC & Programs)	Yes	Yes	Yes	Yes	Yes	Yes
Soft Speech Booster (VAC+)	Yes	Yes	Yes	Yes	No	No
Tinnitus SoundSupport™	Yes	No	Yes	No	Yes	No
Speech Guard E	Yes	Yes	No	No	No	No
YouMatic	Yes (Premium)	Yes (Premium)	Yes (Advanced)	Yes (Advanced)	Yes (Essential)	Yes (Essential)
Free Focus	Yes (Premium)	Yes (Premium)	Yes (Advanced)	Yes (Advanced)	Yes (Essential)	Yes (Essential)
Automatic Directionality	Multi-band Adaptive	Multi-band Adaptive	Multi-band Adaptive	Single Band Adaptive	Single Band Adaptive	Single Band Adaptive
Inium Sense feedback shield	Yes	Yes	Yes	Yes	Yes	Yes
Spatial Sound	Yes (Premium)	No	Yes (Advanced)	No	No	No
Binaural Noise Management	Yes	No	No	No	No	No
Personal Profiles	5	5	3	3	3	3
Noise Management	TriState	TriState	TriState	Modulation-based	Modulation-based	Modulation-based
Special Purpose Programs	Yes	Yes	Yes	Yes	Yes	Yes
ConnectLine/Remote Compatible**	Yes	Yes	Yes	Yes	Yes	Yes
Power Bass & Music Widening (streaming)	Yes	Yes	Yes	Yes	No	No
Fitting Bandwidth (accessed in software)	10 kHz	10 kHz	8 kHz	8 kHz	8 kHz	8 kHz
Fitting Bands	10	10	8	6	6	4
Activity Analyzer (Data Logging)	Yes	Yes	Yes	Yes	Yes	Yes
Life Learning	Yes	Yes	Yes	Yes	No	No
Number of Programs	4	4	4	4	4	4
Fitting Rationales	VAC+, NAL, DSL	VAC+, NAL, DSL	VAC+, NAL, DSL	VAC+, NAL, DSL	NAL, DSL	NAL, DSL
AutoPhone*	Yes (except designRITE)	Yes	Yes (except designRITE)	Yes	Yes (except designRITE)	Yes
T-Coil/FM/DAI (Full-sized BTE & RITE styles)*	Yes	Yes	Yes	Yes	Yes	Yes
FittingLINK Wireless Programming***	Yes	Yes	Yes	Yes	Yes	Yes

*Consult Order Forms for all custom product feature, power, T-coil, style, and wireless options.

**Streamer Pro provides teleloop access and FM compatibility

***Provides wireless programming for all wireless/ConnectLine compatible styles.

Inium Sense Dynamo Matrix

Features	Dynamo SP10	Dynamo SP8	Dynamo SP6	Dynamo SP4
Super Power BTE Style w/13 Battery	Yes	Yes	Yes	Yes
Non-linear fitting formulas	DSE, NAL, DSL	DSE, NAL, DSL	DSE, NAL, DSL	DSE, NAL, DSL
DSE personal profiles	5	3	3	1
Linear fitting formulas	DSElin, NAL-RP	DSElin, NAL-RP	DSElin, NAL-RP	DSElin, NAL-RP
Speech Rescue™	Yes	Yes	Yes	Yes
Speech Guard E	Yes	Yes	No	No
Binaural Noise Management*	Yes	No	No	No
Binaural Synchronization (automatics)	Yes	Yes	No	No
Binaural Coordination (PB operations)	Yes	Yes	Yes	Yes
YouMatic	Premium	Advanced	Essential	No
Transient Noise Management*	Yes	No	No	No
Noise Management	TriState	TriState	Standard	Standard
Inium Sense feedback shield ^{SP}	Yes	Yes	Yes	Yes
Free Focus	Premium	Advanced	Essential	Basic
Automatic Directionality	Multi-band Adaptive	Multi-band Adaptive	Single Band Adaptive	Single Band Adaptive
Power Bass & Music Widening (streaming)*	Yes	No	No	No
Voice Priority i™*	Yes	Yes	Yes	No
Special Purpose Programs	Yes	Yes	Yes	No
Easy RECD™	Yes	Yes	Yes	No
Fitting Bandwidth (accessed in fitting software)	6.5 kHz	6.5 kHz	6.5 kHz	6.5 kHz
Fitting Bands	9	8	6	4
Number of Programs	4	4	4	4
T-Coil, FM, & DAI	Yes	Yes	Yes	Yes
AutoPhone	Yes	Yes	Yes	Yes
FittingLINK Wireless Programming	Yes	Yes	Yes	Yes

*When Speech Rescue™ is enabled, these features are de-activated.

Inium Sense Style Options

IP58




Proven water resistance


All Inium Sense BTE and RITE style instruments take reliability to the next level. Nano-coating of instrument housings and internal components results in a certified rating of IP58.

Oticon hearing instruments are designed to be worn in everyday situations, so users do not need to worry about water or getting wet in a rain shower. The hearing instruments are not intended to be worn while showering or participating in water activities.


Custom styles colors



RITE & BTE styles



Additional colors



RITE



designRITE



miniRITE*



RITE*

miniFit receivers



B0



60



B5



100

(dedicated designRITE receiver)

Power Receiver Molds



100



105

BTE with hook



miniBTE



BTE*



BTE Power



BTE Plus Power

BTE with Corda miniFit thin tube



miniBTE



BTE*



0.9 mm



1.3 mm

Dynamo



BTE 13 Super Power

RITE and BTE miniFit earpieces



Open dome



Bass dome



Power dome



GripTip



LiteTip**



Micro mold**

Custom



IIC



CIC



ITC



Half shell



Full shell

* Available with Tinnitus SoundSupport™ in Alta2 Pro Ti, Nera2 Pro Ti and Ria2 Pro Ti. BTE tinnitus style is only available as BTE85. Power Receiver Mold 105 is not available for the RITE tinnitus styles.

** Also available in VarioTherm material.



www.pro.oticonusa.com

oticon
PEOPLE FIRST