

# Oticon Hearing Aids and DAI/FM Overview

Hearing aids	DAI adapter		FM adapter	FM receiver	Modes of operation	
<p>Automatically adds 2 programs to the hearing aid: - DAI/FM + HA (default) - DAI/FM only</p>	<p>Low Level Input Bodyworn FM solutions e.g. Amigo R5</p>	<p>High Level Input TV, stereo, MP3 player, CD player, computer</p>	<p>Adapter for Amigo receivers and other ear-level receivers</p>	<p>Universal and dedicated Amigo receivers</p>	<p>Default modes of operation when using DAI adapters and FM receivers</p>	
 <p><b>Oticon Alta2/Nera2/Ria2</b> <b>Oticon Alta/Nera/Ria</b> <b>Oticon Dynamo</b> <b>Oticon Sensei</b> <b>Oticon Sensei SP</b> <b>Oticon Safari</b> <b>Oticon Safari SP</b> <b>Oticon Chili</b> <b>Oticon Agil</b> <b>Oticon Acto</b> <b>Oticon Ino</b> <b>Oticon Get**</b> (BTE &amp; RITE styles except designRITE, miniRITE and miniBTE)</p>	<p>AP 900 399-50-520-08 Grey</p> 		<p>FM 9 399-50-590-00 Transparent</p> 	<p>Amigo R2*</p> 	<p><b>Amigo R2 switch positions:</b> ○ or ○○ = HA + FM receiver ● = HA only/FM off</p>	
				<p>Amigo R12 17 colors</p> 		
					<p>Amigo Arc 5 colors</p> 	<p>Via telecoil</p>
				<p>Streamer Pro white 120178 Streamer Pro black 120179</p> 	<p>Amigo R2</p> 	<p><b>Amigo R2 switch positions:</b> ○ or ○○ = HA + FM receiver ● = HA only/FM off</p>
 <p><b>Oticon Opn</b> <b>BTE13 PP (FW 6.0)</b></p> <p><b>Oticon Siya BTE PP &amp; BTE</b></p> <p><b>Oticon Alta2/Nera2/Ria2</b> <b>BTE13 PP</b></p>	<p>AP 1000 142207 Grey</p> 		<p>FM 10 142328 Transparent</p> 	<p>Amigo R2</p> 	<p><b>Amigo R2 switch positions:</b> ○ or ○○ = HA + FM receiver ● = HA only/FM off</p>	
				<p>Amigo R12G2 7 colors</p> 		
					<p>Amigo Arc 5 colors</p> 	<p>Via telecoil</p>
 <p><b>Oticon Sumo family</b></p>	<p>AP 700 399-50-900-08 Grey</p> 		<p>FM 7 399-50-920 Grey</p> 	<p>Amigo R2</p> 	<p><b>Hearing aid positions***:</b> 1 = HA only - M 2 = HA + FM - MT 3 = FM only - T</p> <p><b>Amigo R2 switch positions:</b> ○ or ○○ = HA + FM receiver ● = HA only/FM off</p>	
					<p>Amigo R7</p> 	<p><b>Amigo R7 switch positions:</b> ○ = FM ON ● = FM OFF</p>
	<p>Input levels are set manually</p>				<p>Amigo Arc 5 colors</p> 	<p>Via telecoil</p>




\* On 312 battery style, the LED on Amigo R2 must be set to blink.

\*\* Streamer Pro is not compatible with Oticon Get family.


\*\*\* Other configurations are possible using Genie.

See Amigo  
compatibility guide  
[www.oticon.global/amigo](http://www.oticon.global/amigo)

**oticon**  
PEOPLE FIRST

Cochlear implants	DAI adapter	FM adapter	FM receiver	Modes of operation
 <p>CI processors from various manufacturers</p>		Euro adapter from manufacturer	Amigo R2 	Amigo R2 switch positions: ○ or ○○ = HA + FM receiver ● = HA only/FM off
			Amigo Arc 5 colors 	Via telecoil

## Order information

Amigo Star		
Colors	216-218 MHz	
● C090 Chroma beige	124404	
● C093 Chestnut brown	128692	
● C047 Cool Blue	128695	
● C046 Cool Red	128694	
○ C044 Silver	128693	
● C063 Diamond black	128691	

Amigo R12		Amigo R12G2
Colors	216-218 MHz	216-218 MHz
● Champagne	140-07-100-00	
● C090 Chroma beige	140-07-110-00	162482
● C092 Steel grey	140-07-111-00	188336
● C091 Silver grey	140-07-112-00	188335
● C093 Chestnut brown	140-07-113-00	188332
● C043 Baby pink	140-07-114-00	
● C042 Baby blue	140-07-115-00	
● C047 Cool Blue	140-07-116-00	
● C046 Cool Red	140-07-117-00	
● C034 Purple	140-07-118-00	
○ C044 Silver	140-07-119-00	188334
● C094 Terracotta	102408	188333
● C063 Diamond black	101827	188337
● C048 Emerald green	131234	
● C058 Aquamarine	162068	
● C057 Power pink	162069	
○ C098 Pure white	162070	

Amigo R2		Amigo R7	
Colors	216-218 MHz	216-218 MHz	
● Champagne	140-06-300-00	140-06-600-00	140-06-700-00

DAI Cords		
Cord, Monaural	CD or MP3 Players, computer, etc	383-58-940-09
Cord, Binaural	CD or MP3 Players, computer, etc.	383-58-941-00
DAI Y cord	60 cm (24") white	384-01-066-00
DAI Y Cord	60 cm (24") grey	384-01-063-00
DAI Y Cord	75 cm (30") white	384-01-070-00
DAI Y Cord	75 cm (30") grey	384-01-068-00
DAI Cord	60 cm (24") white	384-01-065-00
DAI Cord	60 cm (24") grey	384-01-062-00
DAI Cord	75 cm (30") white	384-01-069-00
DAI Cord	75 cm (30") grey	384-01-067-00

See Amigo compatibility guide  
[www.oticon.global/amigo](http://www.oticon.global/amigo)