

OTICON | **Opn**

Everything
you need to
stay connected



oticon
PEOPLE FIRST


Make everyday tasks **easier**

Oticon Opn™ wireless accessories are a comprehensive and easy-to-use range of devices developed to improve your listening and communication experiences.

The wireless accessories help you get the most out of your Oticon Opn hearing aids. For example, you can stream the sound from your TV, mobile phone, music player, PC or other devices directly to your hearing aids – like wearing wireless headphones.

Discreet and convenient

For maximum discretion and convenience, you can control your hearing aids with the Oticon ON App or a remote control, connect to the Internet of Things and much more.



Control your hearing aids with Oticon ON App

Oticon ON App lets you choose what you want to hear and adjust the sound. It also offers a low battery notification and more.

Download the free app on the App Store® or Google Play™.



Download on the
App Store



GET IT ON
Google Play

Make clear, safe and **hassle-free** calls

Your Oticon Opn hearing aids can connect to modern mobile phones with Bluetooth™.



Hands-free calls for mobile phones

With ConnectClip, it's easy to make hands-free calls with your iPhone® or Android™ phone.

ConnectClip transmits the caller's voice directly to your hearing aids while the built-in microphone clearly transmits what you say to the caller. ConnectClip transforms your Opn hearing aids into a high-quality wireless headset enabling you to:

- Make and receive hands-free calls in the car or on the go – streaming to both ears
- Wirelessly stream music, videos or other sound from your mobile phone to your hearing aids
- Control the volume and change programs remotely



Made for
 iPod iPhone iPad

Works with
Android™

Oticon Opn is a Made for iPhone® hearing aid

- Connect directly to your iPhone, iPad®, or iPod touch®
- Turn your hearing aids into stereo headphones
- Stream calls, music, videos or other sound wirelessly to your hearing aids
- Use your iPhone to control volume or change programs

Apple, the Apple logo, iPhone, iPad, iPod touch, and Apple Watch are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. Android, Google Play, and the Google Play logo are trademarks of Google Inc. The Android robot is reproduced or modified from work created and shared by Google and used according to terms described in the Creative Commons 3.0 Attribution License.

TV AND STEREO



REMOTE CONTROLS



Get the sound **straight to your hearing aids**

The TV Adapter lets you enjoy your favorite TV program or listen to music from your stereo directly in your hearing aids – just as if you were wearing wireless stereo headphones.

With the TV Adapter you can:

- Clearly hear the TV or listen to music by streaming the sound directly to your hearing aids
- Choose the level you want to hear while leaving the TV at a volume that is comfortable for others
- Use Oticon ON App to control the volume of the streamed sound remotely
- Support Dolby Digital® format
- Range: up to 50 feet



Control your hearing aids **remotely**

If you want to control the volume, change programs or simply turn off your hearing aids with complete discretion, you can use any of the three remote control options: Oticon ON App, ConnectClip or Remote Control.





Hear the **speaker** more easily

Even with the best hearing aids, it can sometimes be difficult to clearly hear what someone is saying – such as in noisy environments or when attending a lecture.

The wireless ConnectClip can help you in difficult listening situations by streaming sound directly to your hearings aids.

Hear speech at a distance or in noise

- Hand ConnectClip to the person you want to hear better
- ConnectClip streams the speaker's voice directly to your hearing aids without affecting your ability to hear other sounds around you
- Use the Oticon ON App to control the volume of the remote microphone
- Range: up to 65 feet

ConnectClip is also ideal for talking over the computer by turning your hearing aids into a wireless headset.


Chat over the computer

- Use ConnectClip to stream the sound directly to your hearing aids when chatting over Skype, Lync, etc.
- Hear the caller's voice directly in your hearing aids while still hearing the sounds around you
- Pair ConnectClip directly to your computer or use Oticon's USB dongle
- Range: up to 33 feet



Up to 33 feet
←-----→





Connect your Opn hearing aids to **the Internet**


The “If This Then That (IFTTT)” network gives you endless ways to wirelessly connect your Oticon Opn hearing aids to other Internet-connected devices. For example, you can:

- Receive notification directly in your hearing aids whenever the doorbell rings
- Have the weather forecast, text messages, breaking news, etc. read out loud in your hearing aids
- Start your coffee maker when you turn on your hearing aids in the morning
- Switch off the lights automatically when you switch off your hearing aids

Find out more about what Oticon Opn™ can do for you. Talk to your hearing care professional today.

Go to **oticon.com/connectivity** for more information on products and compatibility between devices.

Contact your local hearing care professional:



13798 15555-0842/10.17



www.oticon.com

oticon
PEOPLE FIRST