

OTICON | **Opn**

Make clear, **hands-free calls to both ears** from any modern smartphone

Oticon Opn™ enables high-quality streaming with 2.4 GHz Bluetooth Low Energy from any modern smartphone to both hearing aids. Make hands-free calls, stream music, videos or other sounds from iPhone® or Android™ smartphones with the new ConnectClip.



Made for iPod iPhone iPad Works with **Android™**


Apple, the Apple logo, iPhone, iPad, and iPod touch are trademarks of Apple Inc., registered in the U.S. and other countries. Android is a trademarks of Google Inc. The Android robot is reproduced or modified from work created and shared by Google and used according to terms described in the Creative Commons 3.0 Attribution License.

ConnectClip turns Opn hearing aids into a **wireless stereo headset** for use in everyday situations



Remote microphone

Hear what is being said at a distance or in difficult listening environments with ConnectClip's remote microphone functionality.



Computer

Create a wireless headset connection and chat over Skype, Lync, etc. Pair ConnectClip directly to the computer or use Oticon's BTD 800 USB dongle.




Remote control

Use ConnectClip as a remote control to adjust the volume and change programs discreetly with the simple touch of a button.

14018 15555-10180 /10.17

A wide range of Opn wireless accessories with **2.4 GHz high quality streaming**



TV

Oticon's TV Adapter enables sound to be streamed from a TV directly to multiple Opn hearing aids.



Music

Enjoy music in both hearing aids in stereo. Stream directly from iPhone, iPad®, and iPod touch® or use ConnectClip with any Bluetooth-enabled device. Use the TV Adapter for wireless streaming from any music system.



Oticon ON App

Discreetly switch between programs, adjust the sound, receive low battery alerts, use "find my hearing aid" and much more with the Oticon ON App.



Remote control

The Remote Control lets Opn users discreetly regulate volume, switch programs or mute their hearing aids at the touch of a button.

IFTTT Internet connectivity

The If This Then That (IFTTT) network makes life easier and more fun by automatically connecting Opn hearing aids with an endless range of Internet-connected devices used in everyday life.

Learn more about our wireless connectivity accessories at [oticon.com/connectivity](https://www.oticon.com/connectivity)

oticon
PEOPLE FIRST