

OTICON | **Opn**

Opn™ is now rechargeable!

Turn Opn miniRITE into a rechargeable hearing aid



The world's only rechargeable hearing aid with an open sound experience

- Open sound experience: delivered by OpenSound Navigator™ and Spatial Sound™ LX
- BrainHearing™ benefits: 20% less listening effort, 20% more capacity to remember, 30% better speech understanding*
- 96% of users hear better or much better with Opn**


The world's first rechargeable hearing aid with 2.4 GHz direct streaming

- World's first dual radio communication with TwinLink™ technology: 2.4 GHz and NFMI
- Made for iPhone® hearing aid
- World's first Internet-connected hearing aid

* Le Goff et al. 2016

**700 Opn users voluntarily mailed response cards to share their personal perceptions of the open sound experience. The users (average age 61-70 years) had worn Opn for 12 weeks or less.

Apple, the Apple logo, iPhone, iPad, iPod touch, and Apple Watch are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. Android, Google Play, and the Google Play logo are trademarks of Google Inc.



oticon
PEOPLE FIRST

A Powerful Combination

Oticon Opn miniRITE rechargeable option is easy to use and very convenient. No need to replace batteries – the hearing aids are charged overnight for a full day's use.

Rechargeable hearing aids

- **Convenient:** no need to change batteries every few days.
- **Reassuring:** fully charges overnight to ensure sufficient power throughout the day.
- **Upgradable:** turn any client's Opn miniRITE into a rechargeable device by exchanging the battery door.



Charger

- **Easy to use:** simply place the hearing aids in the charger overnight.

Rechargeable batteries

- **Practical:** saves time buying and disposing of hundreds of batteries a year.
- **Maximum flexibility:** the 312 ZPower batteries can be exchanged with disposable batteries if your client forgets to charge them overnight.
- **Environmentally-friendly:** one pair of ZPower batteries can save around 150-200 disposable batteries a year.

Estimated Battery Life

No streaming	18-19.5 hours
90 minutes of iPhone Streaming	16-17 hours
4 hours of TV +1 hour of iPhone streaming	15-16 hours

Battery Life will vary depending upon hearing loss and hearing aid use.