

Oticon Opn™ Style Options

miniRITE, miniRITE-T, BTE13 PP



Oticon Opn miniRITE – small, discreet, with single push button for volume and program changes.

Oticon Opn miniRITE-T – discreet style featuring a telecoil and convenient double push button for easy volume and program control.

Oticon Opn BTE 13 PP – compact power style with tactile double push button and LED light for monitoring hearing aid status.


All styles offer 2.4 GHz wireless technology and made for iPhone® functionality.*

miniFit Receivers and miniFit Corda Thin Tube options

Select between three receivers with different output performance, with lengths 0-5 for non power speakers and 1-5 for power.



Corda miniFit thin tube



miniRITE and Corda Thin Tube Earpieces

miniFit earpieces include standard domes, MicroMold, and LiteTip. Custom earpieces are available in acrylic, soft silicone or VarioTherm® material.



Open dome Bass dome, single vent (0.8 mm) Bass dome, double vent (1.4 mm) Power dome Grip Tip No vent Grip Tip Vent MicroMold LiteTip MicroMold, VarioTherm® LiteTip, VarioTherm®

Power Receiver Molds



IP68

Proven IP68 water and dust resistant hearing aids.



12613 15555-10020/3.17



www.oticon.com



OTICON | **Opn**

Product comparison

	Oticon Opn 1	Oticon Opn 2	Oticon Opn 3
Speech Understanding			
OpenSound Navigator™	Level 1	Level 2	Level 3
Balancing power effect	100%	50%	50%
Max. noise removal	9 dB	5 dB	3 dB
Speech Guard™ LX	Level 1	Level 2	Level 3
Spatial Sound™ LX	4 estimators	2 estimators	2 estimators
Soft Speech Booster LX	•	•	•
Sound Quality			
Speech Rescue	•	•	•
Clear Dynamics	•	•	-
Binaural Noise Management	•	•	-
Fitting Bandwidth (accessed in software)	10 kHz	8 kHz	8 kHz
Processing Channels	64	48	48
Bass Boost (streaming)	•	•	•
Listening Comfort			
Transient Noise Management	4 configurations	On/Off	On/Off
Feedback shield LX	•	•	•
Wind Noise Management	•	•	•
Binaural Coordination	•	•	•
Personalization & Optimizing Fitting			
YouMatic™ LX	3 configurations	2 configurations	1 configuration
Fitting Bands	16	14	12
Multiple Directionality Options	•	•	•
Adaptation Management	•	•	•
Oticon Firmware Updater	•	•	•
Tinnitus SoundSupport™	•	•	•
Speech Rescue	•	•	•
Fitting Formulas	VAC+, NAL, DSL	VAC+, NAL, DSL	VAC+, NAL, DSL
Acoustic Notifications	•	•	•
Connecting to the World			
Stereo streaming (2.4 GHz)	•	•	•
Made for iPhone®	•	•	•
Oticon ON App	•	•	•
ConnectClip (coming in 2017)	•	•	•
Remote Control 3.0	•	•	•
TV Adapter 3.0	•	•	•
AutoPhone	•	•	•
Special Needs			
Tinnitus SoundSupport™	•	•	•

*Oticon Opn™ is compatible with iPhone 7 Plus, iPhone 7, iPhone SE, iPhone 6s Plus, iPhone 6s, iPhone 6, iPhone 5s, iPhone 5c, iPhone 5, 9.7-inch iPad Pro, 12.9-inch iPad Pro, iPad Air 2, iPad Air (4th generation), iPad mini 4, iPad mini 3, iPad mini 2, iPad mini, and iPod touch 5th and 6th generation. Devices must be running iOS 9.3 or later. Please visit www.oticon.com for more details.