

Ship To Information

Customer Number:

(Please complete all information including name & phone number)

Phone #:() _____ Purchase Order #: _____

Company Name: _____

Address: _____

City: _____ State: _____ Zip: _____

Fitter's Information

Today's Date: _____ Fitting Date: _____

Fitter's Name: _____

Fitter's E-mail: _____

Patient Information

First Name: Middle Initial: Age:

Last Name:

Last 4 Digits of Social Security #:

Bill To Information

Bill To Number:

Alta2 Pro

	Alta2 Pro	Battery	Power (90)	Ultra Power (100)	Omni Only	PB*†	VC	Auto Phone (not available with Tcoil)	T-coil	Auto-T	Remote Control	Streamer Pro 1.3A
		L R	L R	L R	L R	L R	L R	L R	L R	L R	Black	Black White
IIC**		10A <input type="checkbox"/> <input type="checkbox"/>	NA	NA	STD	NA	NA	NA	NA	NA	NA	NA
CIC		10A <input type="checkbox"/> <input type="checkbox"/>	<input type="checkbox"/> <input type="checkbox"/>	NA	STD	<input type="checkbox"/> <input type="checkbox"/>	NA	NA	NA	NA	<input type="checkbox"/>	<input type="checkbox"/> <input type="checkbox"/>
Canal		312 <input type="checkbox"/> <input type="checkbox"/>	<input type="checkbox"/> <input type="checkbox"/>	NA	<input type="checkbox"/> <input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> <input type="checkbox"/>					
Half Shell		312 <input type="checkbox"/> <input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> <input type="checkbox"/>								
		13 <input type="checkbox"/> <input type="checkbox"/>	<input type="checkbox"/> <input type="checkbox"/>	<input type="checkbox"/> <input type="checkbox"/>	NA	<input type="checkbox"/> <input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> <input type="checkbox"/>				
Full Shell		312 <input type="checkbox"/> <input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> <input type="checkbox"/>								
		13 <input type="checkbox"/> <input type="checkbox"/>	<input type="checkbox"/> <input type="checkbox"/>	<input type="checkbox"/> <input type="checkbox"/>	NA	<input type="checkbox"/> <input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> <input type="checkbox"/>				

All models default to Directional Microphones unless Omni is selected.

*Push button needed for T-coil access. Push button is default on custom instruments unless otherwise indicated.

**non-wireless

†ConnectLine™ Options

- TV Adaptor
- Phone Adaptor
- ConnectLine Mic

‡ Custom Instruments offer only 1 control per instrument.

Select one option for each instrument if desired.

- Left VC OR PB
Right VC OR PB

More Options

Color Options

- | | | |
|--|---|---|
| IIC Faceplate | IIC Shell | Alta Pro |
| <input type="checkbox"/> Black (default) | <input type="checkbox"/> Red/Blue (default) | <input type="checkbox"/> Pink (default) |
| <input type="checkbox"/> Pink | <input type="checkbox"/> Pink | <input type="checkbox"/> Pink (clear shell) |
| <input type="checkbox"/> Tan | <input type="checkbox"/> Clear | <input type="checkbox"/> Tan |
| <input type="checkbox"/> Med. Brown | | <input type="checkbox"/> Tan (clear shell) |
| <input type="checkbox"/> Dark Brown | | <input type="checkbox"/> Medium Brown (clear shell) |
| | | <input type="checkbox"/> Dark Brown (clear shell) |

Audiometric Information:

Hz	250	500	1K	2K	3K	4K	8K
AC Right							
AC Left							

Impression Information:

Type: Open jaw Closed jaw Return Impressions

Canal Length: (Factory STD unless otherwise marked)

- Medium Long As marked on impressions

Venting: Factory Selected (largest possible)

- Style:** Small Vent Large Vent SAV Vent No Vent
Vent Diameter: Medium Vent Max Vent Pressure Vent

Other Options

- Raised VC Canal Lock
- Removal String (Standard on IIC/CIC) Raised PB
- Clothing Loop Removal notches
- Include Spanish Instructions

Special Instructions

Venting

All vents are calculated in Genie, unless other wise indicated on the order form. For a specific vent type please note it in special instructions.



Collection Vent
Cone shape channel, round opening faceplate, "D" opening on canal tip. Most common for ITC, HS, and FS where pressure vents not necessary.



Cast Vent
Straight channel through device. Round opening on faceplate and canal tip. Used when pressure vents are indicated.



Oval Vent
Oval openings used for smaller instruments IIC and CIC.



Free Vent
Free form openings, takes advantage of space to ensure largest vent possible.

Accessories



T-Cap
The nano-coated, wax-repellent T-Cap filter is used for all miniature custom styles.



O-Cap
The new nano-coated O-Cap microphone protection filter is designed to reject wax and fits all instruments using battery size 312 and 13.



ProWax
Providing advanced receiver port protection, the ProWax filter fits all styles in the new Oticon custom program.



T Cap (mic cover)
Gray compact/Black tools

Part Number (1 pc):
123328 - Beige
123330 - Light Brown
123331 - Med Brown
123329 - Dark Brown
123327 - Black



O Cap (mic cover)
Gray compact/Black tools

Part Number (1 pc):
128003 - Beige
128004 - Light Brown
128005 - Med Brown
128006 - Dark Brown



ProWax
Gray "turtle"/Lt Gray tools

Part Number (1 pc):
123367



Custom Option for Special Needs

Instrument Removal Options	Large Ball 	Reinforced Removal Line 	Door Hook 
Tactile Options	Raised Push Button 	VW Large, Small or Raised 	Battery Indicator 
Clothing Loops	Loops: Metal Small and Large 	Plastic: Small and Large 	Plastic Combo: Loop+removal string 
Retention Options	Canal lock 	Helix Lock 	Spiral Retention Ring 

Please include your request on the comment line of the order form. Options available space permitting and may not be available on all models.

*Push button serves as program change; T-coil access and volume change.

†All models default to Directional Microphones unless Omni is selected.

‡Technical information re Zcc coupler:
Power Receiver: OSPL90 PK 121, PK Gain: 54, Ultra Power Receiver: OSPL90 PK 127, PK Gain: 62