

PRODUCT INFORMATION REMOTE CONTROL 2.0

*The **Remote Control 2.0** allows users to operate Oticon hearing instruments discreetly without any other devices needed.*

The Remote Control is a convenient solution for all wireless instruments. It is an excellent option for users whose instruments do not have controls.

Easy to operate with full discretion

The Remote Control offers a straightforward and discreet way to operate the hearing instruments.

With the Remote Control, the hearing instrument user can perform the most common adjustments without drawing attention to the hearing instrument.

Users can adjust volume up/down and cycle through the hearing instrument programs.

Fully discreet, the Remote Control can operate the hearing instruments from a distance of up to 1 meter; a direct line of sight between the Remote and the hearing instruments is not required.

Easy to operate

The Remote Control buttons are clearly visible and easy to operate.

Each time a button is pressed, a discreet clicking sound is heard and felt by the user, while a visual indication on the Remote Control confirms the action.



PRODUCT FEATURES



Technical data Remote Control 2.0

Weight (incl. batteries)	32.7 g
Battery Type	2 x AAAA
Battery Life Time at Normal Use	1 year
Typical Transmission Range	1 meter
Indicator	Green LED blinking when Remote Control is operated or battery is low
Operating temperature	10 - 45 °C
Operating humidity	Relative humidity < 85%
Storage temperature	-20 - 70 °C